

4Closure™

SURGICAL FASCIA CLOSURE SYSTEM




FOR FAST, EASY AND SAFE CLOSURE OF TROCAR WOUNDS

The new disposable, patent pending 4Closure™ system includes a suture delivery needle, a needle guide and it features a new operating method for closure of the fascia in laparoscopic procedures that use large operating ports or trocars.

The device can be applied through the operating cannula or can be inserted directly into the wound.

The 4Closure™ kit is intended to prevent post operative hernia defects and is also used to help control acute bleeding from the abdominal wall. Both ends of the 4Closure™ guide rod are usable for closing the wound depending on the patient's weight.

Available sizes: 10-11mm, 12mm, 15mm and 23mm



Suture Passer Instrument

OPEN



INLINE SLIDING MECHANISM



CLOSED

BENEFITS:

- Compatible with existing trocars
- No loss of pneumoperitoneum
- Fast closure of the fascia
- Rapid learning curve
- Suitable for variety of laparoscopic applications
- No risk of mechanical failure during surgery
- Crucial for Laparoscopic Bariatric Surgery
- Major decrease in post-surgical pain in trocar site by eliminating the bundling of tissue, muscle and fascia

For more information please contact:

Allez Medical Applications, Inc.
2141 South Standard Ave.,
Santa Ana, CA 92707

Phone: 714.641.2098
Fax: 714.557.5361
E-mail: info@allezmedical.com
Web: www.allezmedical.com

